



Maja Krstic Ristivojevic

Maja Krstic Ristivojevic

Microplastics

18. Septembar 2025. Novi Sad



Introduction


EXploration of **P**ETase side **A**ctivity of digestive enzymes of human gastrointestinal tract acting on micro- and nanoplastics: mode of action and products **Charac**Terization

Acronym: **XPACT**

Participating Scientific and Research Organizations (SROs):

- * University of Belgrade, Faculty of Chemistry (UBFC)
- * Institute of Chemistry, Technology and Metallurgy (ICTM)

Principal Investigator (PI): Tanja Cirkovic Velickovic




Good neighbours
creating
common future



Introduction

- Microplastic size distribution: mesoplastics (5 mm - 25 mm), microplastics (1 µm - 5 mm), and nanoplastics (less than 1 µm).
- Polyethylene terephthalate - Terephthalic acid (TPA) + Ethylene glycol (EG)




- Enzymatic degradation of PET
- Hydrolases - esterases, lipases

*Good neighbours
creating
common future*



Introduction - FTIR Spectroscopy (Fourier Transform Infrared Spectroscopy)



Purpose:

- Identify chemical bonds and functional groups in molecules
- Analyze solids, liquids, or gases

Principle:

- Molecules absorb infrared (IR) light at specific wavelengths
- Absorption causes vibrational transitions → “molecular fingerprint”

ATR-FTIR (Attenuated Total Reflectance)

- Specialized FTIR technique
- IR light reflects internally on the crystal in contact with the sample
- Measures the surface of the sample directly without complex preparation

How it Works:

- IR light passes through the ATR crystal into the sample
- Detector measures absorbed wavelengths
- Fourier Transform converts raw data to IR spectrum (absorbance vs. wavenumber)

Key Features:

- Peaks correspond to functional groups (e.g., O–H, C=O, C–H)
- Non-destructive, fast, minimal sample preparation
- Widely used in chemistry, materials, forensic science, and biochemistry

Good neighbours
creating
common future



Microplastics




Figure 1. Fluorescence images of PET MPs with (Ctrl+) or without (Ctrl-) BSA-AF488 hard corona acquired with 10 x (A) and 4 x magnification (B), and after staining with Nile red (C) following treatment with three clean-up protocols: I (10% SDS + 15% H₂O₂); II (2 x 30% H₂O₂); III (15% H₂O₂ + 10% KOH). Images were acquired with 10x magnification after excitation at 340–380 nm (UV filter for BSA-AF488) and/or 527.5–552.5 nm (G – green filter for Nile Red).





Figure 2. Fluorescence intensity (%) (A), and fluorescent image of the gel after reducing SDS-PAGE (B) of PET MPs with (+) or without (-) BSA-AF488 hard corona following treatment with three clean-up protocols: I (10% SDS + 15% H₂O₂); II (2 × 30% H₂O₂); III (15% H₂O₂ + 10% KOH). Ctrl+ - Positive control (PET MPs with BSA-AF488 hard corona without treatment). The fluorescence intensity was normalized to the mass of the weighed PET MP particles. 30 µL of each sample was loaded into each well.

Microplastics






Figure 3. ATR-FTIR spectra of PET MPs with BSA-AF488 hard corona compared to untreated PET MPs.



A



B




Figure 4. ATR-FTIR spectra (A) and overlaid spectra (B) of PET MPs incubated without BSA-AF488 and treated with different clean-up protocols: I- (10% SDS + 15% H_2O_2), II- ($2 \times 15\%$ H_2O_2), and III- (15% H_2O_2 + 10% KOH), compared to untreated PET MPs (Ctrl-).

Good neighbours
creating
common future

Microplastics





Figure 5. Overlaid spectra (left) and the difference spectra (right) of PET MPs incubated: 24h with chymotrypsin 25 IU/ml (A), 24h with lipase 930 IU/mL (B) and 72h with lipase 75.6 IU/mL in simulated intestinal fluid and corresponding controls (C) (PET MPs incubated without enzymes under the same conditions). The clean-up protocol III (H_2O_2 + KOH) was used for removing hard and soft coronas.

Good neighbours
creating
common future



Microplastics



ENZYME	ACTIVITY	
Lipase from porcine pancreas, cat.no. L3126-100G, pcode 1003455012, Source SLCN3698	28 U/mg	
Pepsin (XPACT), cat.no. P7012-1G, Source SLCR2106	1418 ± 43 U/mg	
	LIPASE ACTIVITY	TRYPSIN ACTIVITY
Pancreatin (IMPTOX), cat.no. P1750-100G, Pcode 1003351803, source SLCK2763 (solubility in water 7 mg/ml)	25.62 U/mg	1.13 ± 0.07 U/mg
Pancreatin, cat.no. P7545-100G, lot 0000344645, 0000324607 (solubility in water 20 mg/ml)	96.8 U/mg	6.99 ± 0.58 U/mg

Figure 6. SDS-PAGE of digestive enzymes used in the experiments on 14% polyacrylamide gel. Mw – molecular weight marker; T – trypsin; C – chymotrypsin; L-lipase.



Thank you!

• **Тамара Д. Васовић**

Универзитет у Београду - Хемијски факултет (Београд)

• **Стефан Г. Ивановић**

Научна установа Институт за хемију, технологију и
металургију (Београд)

• **Весна Б. Јовановић**

Универзитет у Београду - Хемијски факултет (Београд)

• **Маја В. Крстић Ристивојевић**

Универзитет у Београду - Хемијски факултет (Београд)

• **Тамара Н. Лујић**

Универзитет у Београду - Хемијски факултет (Београд)

• **Тамара Р. Мутић**

Универзитет у Београду - Хемијски факултет (Београд)

• **Маја Младеновић Стоканић**

Универзитет у Београду - Хемијски факултет (Београд)

• **Ана С. Симовић**

Универзитет у Београду - Хемијски факултет (Београд)

• **Драгана Ј. Станић-Вучинић**

Универзитет у Београду - Хемијски факултет (Београд)

• **Тања Д. Ћирковић Величковић**

Универзитет у Београду - Хемијски факултет (Београд)

*Good neighbours
creating
common future*